

Optimal Exposure Minimal Trauma



Unique curved rack follows anatomic curves of the body

Opens with a gentle radius and minimizes trauma to the ribcage




Genesee
BIOMEDICAL, INC.

Design beyond standard.

Adult Sternal Retractor with Interchangeable Atrial and Aortic Rakes

Multipurpose Retractor for Valve Repair or Replacement and On or Off-Pump CABG



Retractor with Valve Exposure Rakes



Retractor with Minimally Invasive Blades

Innovative Design

- Unique rack design opens with a gently curved radius, optimizing exposure with minimal trauma to the chest and ribcage.
- Anatomic curved shape distributes opening pressure differentially on the sternum.
- Pivoting blades equalize pressure points across the sternum.

Versatile System

- Unique ball clamps, H bar mounting rods and pillar clamps provide fast set up.
- Interchangeable rakes provide excellent exposure and ease of use.
- Changeable blades allow adjustment for individual patient anatomy.

Minimally Invasive Approaches

- Curved rack with hinged arms hug patient's body unobtrusively.
- Deep or standard minimally invasive blades optional.
- Optimal platform for off-pump procedures.

Ordering Information

CATALOG NUMBER	PRODUCT DESCRIPTION
SRCL-AR	Sternal retractor with latch system and atrial and aortic rakes
SRCL-1	Sternal retractor with latch system
ST-ST5	Standard sterilization tray
AAR-1	Aortic and atrial rake attachments

Genesee BioMedical reserves the right to change or modify specifications without notice. U.S. Patents #5,772,583 and #6,206,828.
© 2011 Genesee BioMedical, Inc. All Rights Reserved. Printed in USA