




## CT DIVISION



## TeDan Surgical Innovations

Introduction forceps clip onto  
Atrial Lift Blades and introduce  
them into the thoracic cavity  
through the thoracotomy



CT-0242  
**Blade Introduction Forceps**



CT-0240  
**Blade Holder with 5mm Shaft**  
Secures lift blade via intercostal puncture



CT-0210 through CT-0214  
**Atrial Lift Blades**

- 35mm width
- Lengths available in 10mm increments from 40-60mm




CT-0220 through CT-0224  
**Atrial Lift Blades**

- 45mm width
- Lengths available in 10mm increments from 50-70mm

# Atrial Lift System **CT**

DESIGNED TO MEET THE NEEDS OF THORACIC ACCESS



**ML-0039**

**Phantom ML Articulating Accessory Arm**  
with Quick Connector

**ML-0043**

**Phantom ML Instrument Holder**  
with Quick Connect for Articulating Arm



## **CT-1010** ATRIAL LIFT SYSTEM

CAT NO.	DESCRIPTION
CT-0210	Atrial Lift Blade 35W x 40L
CT-0212	Atrial Lift Blade 35W x 50L
CT-0214	Atrial Lift Blade 35W x 60L
CT-0220	Atrial Lift Blade 45W x 50L
CT-0222	Atrial Lift Blade 45W x 60L
CT-0224	Atrial Lift Blade 45W x 70L
CT-0240	Blade Holder with 5mm Shaft
CT-0710	Atrial Lift Sterilization Tray
CT-0242	Blade Introduction Forceps
ML-0039	Articulating Accessory Arm Straight Column with Quick Connector
ML-0041	Table Clamp for Articulating Arm
ML-0043	Phantom ML Instrument Holder with Quick Connect for Articulating Arm

## CONTACT US

TeDan Surgical Innovations, LLC  
 12615 W. Airport Blvd, Suite 200 • Sugar Land, TX 77478  
 Ph: 713-726-0886 • Fax: 713-726-0846 • Toll Free: 877-726-0886  
 www.tedansurgical.com • sales@tedansurgical.com

© TSI 2012

