


Adult Curved Atrial Rakes

See beyond the curve
with new curved atrial
rakes.



Expanded Exposure

- Curved Atrial Rakes for Genesee Biomedical Adult Sternal Retractor
- Less obtrusive, provides optimal exposure curving out of the field of visibility
- Interchangeable rakes combined with ball clamps and H bar produce excellent tissue exposure
- Unique shapes of ball and pillar clamps reduce risk of suture snagging



Genesee
BIOMEDICAL, INC.

Design beyond standard.



Curved Rake, Ball Clamp and H Mounting Bar



Ball Clamps and Mounting Bars

Ordering Information

Catalog Number

Product Description

Adult Curved Atrial Rakes – Standard Handles (6.25" Long)

AR-SLRC	Small Curved Left Atrial Rake (0.75" wide by 1.4" long by 1" deep)
AR-MLRC	Medium curved Left Atrial Rake (1.25" wide by 1.4" long by 1.1" deep)
AR-LLRC	Large Curved Left Atrial Rake (1.5" wide by 2.4" long by 1.75" deep)
AR-LRRC	Large Curved Right Atrial Rake (1.75" wide by 1.6" long by 0.7" deep)
AR-LARC	Long Curved Atrial Rake (0.75" wide by 2.5" long by 2" deep)

Adult Curved Atrial Rakes-Long Handles (8.25" Long)

AR-SLRLC	Small Left Atrial Rake (0.75" wide by 1.4" long by 1" deep)
AR-MLRLC	Medium Left Atrial Rake (1.25" wide by 1.4" long by 1.1" deep)
AR-LLRLC	Large Left Atrial Rake (1.5" wide by 2.4" long by 1.75" deep)
AR-LRRLC	Large Left Atrial Rake (1.75" wide by 1.6" long by 0.7" deep)
AR-LARLC	Long Atrial Rake (0.75" wide by 2.5" long by 2" deep)

Accessories for Curved Atrial Rakes

AR-HW	Large Hex Wrench
AR-HWG	Large Hex Wrench - Green Handle
AR-CC	Adult Ball Clamp for Curved and Straight Rakes
AR-CCT	Adult Toggle Ball Clamp for Curved and Straight Rakes
AR-PC	Adult Pillar Clamp
AR-PCT	Adult Pillar Clamp with Toggle
AR-HR	"H" Mounting Bar
AR-MB	Straight Mounting Bar



Genesee
BIOMEDICAL, INC.

Design beyond standard.

Genesee BioMedical reserves the right to change or modify specifications without notice.
© 2013 Genesee BioMedical, Inc. All Rights Reserved. Printed in USA

www.geneseebiomedical.com