



**CT ULTRAVISION™**

# MITRAL VALVE EXPOSURE

## SYSTEMS

PROVIDING OPTIMAL VISUALIZATION  
FOR MINIMALLY INVASIVE MITRAL VALVE  
PROCEDURES



**TeDan Surgical  
Innovations**

# CT™ MITRAL VALVE EXPOSURE


## ATRIAL LIFT SYSTEM: CT-1010

Retraction of Anterior Atrial Wall



CT-0242  
**Blade Introduction Forceps**

- Introduction Forceps used to introduce Atrial Lift Blades through the thoracotomy incision



CT-0240  
**Blade Holder with 5mm Shaft**

- Secures lift blade via intercostal puncture



CT-0210 - CT-0214  
**Atrial Lift Blades**

- 35mm width
- Lengths available in 10mm increments from 40 to 60mm




CT-0220 - CT-0224  
**Atrial Lift Blades**

- 45mm width
- Lengths available in 10mm increments from 50 to 70mm

## AD ATRIAL RAKE RETRACTOR: CT-0245




Retraction of Anterior Atrial Wall



CT-0245  
**AD Atrial Rake Retractor**

- Introduced through the thoracotomy incision

## POSTERIOR ATRIAL RETRACTOR (PAR): CT-0305



CT-0305

### Posterior Atrial Retractor

- Integrates standard sump line to maintain clear field
- Connects anywhere on retractor rack



LS-0451

### Phantom™ T-Handle

- Used to tighten flex arm on CT-0305
- Same handle used for pivoting mechanism on CT Ultravision™ Retractor

## TABLE MOUNTED ARMS AND TABLE CLAMP



ML-0043

### Phantom ML™ Instrument Holder

ML-0039

### Phantom ML™ Articulating Accessory Arm w/ Quick Connect


- Alternative option to ML-0038 Arm
- Provides multiple attachment options



ML-0038

### Phantom ML™ Articulating Accessory Arm w/ Instrument Holder

- Can be used with CT-1010 or CT-0245



ML-0041

### Phantom ML™ Table Clamp

## TeDan Surgical Innovations, LLC

Ph: 713-726-0886 Fax: 713-726-0846 Toll Free: 877-726-0886

12615 West Airport Blvd., Suite #200, Sugar Land, Texas 77478

Email: [sales@tedansurgical.com](mailto:sales@tedansurgical.com)

[www.tedansurgical.com](http://www.tedansurgical.com)

## MITRAL VALVE EXPOSURE PRODUCT OFFERING

CAT. NO.	DESCRIPTION
CT-0210	ATRIAL LIFT BLADE 35W x 40L
CT-0212	ATRIAL LIFT BLADE 35W x 50L
CT-0214	ATRIAL LIFT BLADE 35W x 60L
CT-0220	ATRIAL LIFT BLADE 45W x 50L
CT-0222	ATRIAL LIFT BLADE 45W x 60L
CT-0224	ATRIAL LIFT BLADE 45W x 70L
CT-0240	BLADE HOLDER w/5MM SHAFT TO SECURE LIFT BLADE VIA INTERCOSTAL PUNCTURE
CT-0242	BLADE INTRODUCTION FORCEPS FOR INSERTION OF LIFT BLADES THROUGH THORACOTOMY
CT-0245	AD ATRIAL RAKE RETRACTOR
CT-0305	POSTERIOR ATRIAL RETRACTOR
CT-0701	ULTRAVISION CT ATRIAL LIFT STERILIZATION TRAY
LS-0451	PHANTOM T-HANDLE FOR RETRACTOR WITH PIVOTING ARMS
ML-0038	ARTICULATING ACCESSORY ARM STRAIGHT COLUMN W/INSTRUMENT HOLDER
ML-0039	ARTICULATING ACCESSORY ARM STRAIGHT COLUMN W/QUICK CONNECTOR
ML-0041	TABLE CLAMP FOR ARTICULATING ARM
ML-0043	PHANTOM ML INSTRUMENT HOLDER W/QUICK CONNECT F/ARTICULATING ARM



CT-0701  
Atrial Lift Sterilization Tray