

 **CT LIFT SYSTEM™**
REINVENTING IMA ACCESS




*For Additional Information
Contact PEAK MEDICAL
info@peakmdi.com
630-917-8234*

TSI TeDan Surgical
Innovations

CT LIFT SYSTEM

IMA Exposure

TSI is reinventing IMA access for standard and minimally invasive procedures by improving on proven technology with a modern design that provides optimal exposure and more positioning versatility. Cut the cable – the sleek design eliminates use of braided steel cables which simplifies cleaning and reduces required maintenance.



Sternal Setup Using Dual Hook

Open IMA Exposure

CT-0052

Dual Hook

- Time tested design to engage the sternum
- Improved connection to lifting mechanism

MI Thoracotomy Setup Using IMA Blade

Minimally Invasive IMA Approach

CT-0050, CT-0300 -CT-0302

Blade Hook and IMA Blades

- Wide and narrow IMA blades attach to existing CT-0100 Ultravision™ retractor accommodating diverse patient anatomies
- Pivot function of retractor improves access and visibility

Versatility Made Simple

Adaptive setup options ensure patient-specific compatibility

Unique adjustment mechanism of extension arm enables various angles and post lengths, providing exact positioning.



Improving System Assembly

Creating effortless and stable IMA exposure

CT Lift System components simplify the setup process through a unique table locking post and one-step cam mechanisms.



Cam mechanism provides precise positioning along the posts and creates a rock-solid lock in one turn of the handle

Universal hook holder optimizes attachment of preferred IMA approach to the lift mechanism



- Table post attaches directly to OR side rail over the drape
- Tightening handle stows when not in use
- No stainless steel braided cables

CT-1705 CT LIFT SYSTEM


Reinventing IMA Access:

- For standard and minimally invasive approaches
- Rapid setup
- Easy to clean
- Lifetime warranty

*For Additional
Information
Contact PEAK MEDICAL
info@peakmdi.com
630-917-8234*

Set includes:

Cat. No.	Description
CT-0050	IMA Blade Hook
CT-0052	Open IMA Dual Hook
CT-0300	IMA Blade, Wide
CT-0302	IMA Blade, Narrow
CT-0705	Sterilization Tray, CT-1705
HS-0016	Table Clamp, 450 mm, Purple
HS-0019	Extension Arm
HS-0100	Lift Assembly, Purple



TeDan Surgical Innovations, LLC

Ph: 713-726-0886 Fax: 713-726-0846 Toll Free: 877-726-0886

12615 West Airport Blvd., Suite #200, Sugar Land, Texas 77478

Email: info@tedansurgical.com

www.tedansurgical.com

LIT-0089, REV 3 ©TSI 2017