



**TeDan Surgical
Innovations**

CT RETRACTOR BLADES

DESIGNED TO MEET THE DIVERSE DEMANDS
OF THORACIC ACCESS

THORACOTOMY BLADES

Concave and convex blades are shaped to conform to the anatomy the ribs.

- Reduces potential for trauma.
- Increases stability in the incision.

GT window blades are now available.

Ultravision CT Retractor Blades are designed in multiple lengths, widths and styles to accommodate diverse patient anatomies. Dual light ports enable the use of TeDan Surgical's LED lighting system.

All aluminum blades are radiolucent and ideal for procedures using fluoroscopy.

*For Additional Information
Contact PEAK MEDICAL
info@peakmdi.com
630-917-8234*



STERNAL BLADES

- Flat blades designed for Mini-Sternotomy access.



PORT BLADES

- "Micro" (30mm) width is ideal for port sized incisions and robotic surgeries.





**TeDan Surgical
Innovations**

CT RETRACTOR BLADES

ORDERING INFORMATION



Stainless steel thoracotomy blades available

THORACOTOMY BLADES

CAT. NO.	DESCRIPTION
CT-0110	XSMALL BLADE CURVED, 45X35MM, GREEN
CT-0111	XSMALL BLADE UPPER RIB (REVERSE CURVE), 45X35MM, GREEN
CT-0112	XSMALL BLADE CURVED, 45X35MM W/ WINDOW, GREEN
CT-0113	XSMALL BLADE UPPER RIB (REVERSE CURVE), 45X35MM W/ WINDOW, GREEN
CT-0114	XSMALL BLADE CURVED, 45X35MM W/LIGHT CHANNEL, GREEN
CT-0115	XSMALL BLADE UPPER RIB (REVERSE CURVE), 45X35MM W/LIGHT CHANNEL, GREEN
CT-0120	SMALL BLADE CURVED, 45X50MM, RED
CT-0121	SMALL BLADE UPPER RIB (REVERSE CURVE), 45X50MM, RED
CT-0122	SMALL BLADE CURVED, 45X50MM W/ WINDOW, RED
CT-0123	SMALL BLADE UPPER RIB (REVERSE CURVE), 45X50MM W/ WINDOW, RED
CT-0124	SMALL BLADE CURVED, 45X50MM W/LIGHT CHANNEL, RED
CT-0125	SMALL BLADE UPPER RIB (REVERSE CURVE), 45X50MM W/LIGHT CHANNEL, RED
CT-0130	MEDIUM BLADE CURVED, 45X70MM, BLUE
CT-0131	MEDIUM BLADE UPPER RIB (REVERSE CURVE), 45X70MM, BLUE
CT-0132	MEDIUM BLADE CURVED, 45X70MM W/ WINDOW, BLUE
CT-0133	MEDIUM BLADE UPPER RIB (REVERSE CURVE), 45X70MM W/ WINDOW, BLUE
CT-0134	MEDIUM BLADE CURVED, 45X70MM W/LIGHT CHANNEL, BLUE
CT-0135	MEDIUM BLADE UPPER RIB (REVERSE CURVE), 45X70MM W/LIGHT CHANNEL, BLUE

STERNAL BLADES

CAT. NO.	DESCRIPTION
CT-0140	FLAT BLADE - STERNAL, 50X60MM
CT-0141	FLAT BLADE - STERNAL, 50X60MM W/ WINDOW
CT-0142	FLAT BLADE - STERNAL, SMALL PATIENT, 50X40MM
CT-0143	FLAT BLADE - STERNAL, SMALL PATIENT, 50X40MM W/ WINDOW
CT-0144	FLAT BLADE - STERNAL, 50X60MM W/LIGHT CHANNEL, ORANGE
CT-0145	FLAT BLADE - STERNAL, SMALL PATIENT, 50X40MM W/LIGHT CHANNEL, ORANGE
CT-0146	FLAT BLADE - STERNAL, SHALLOW, 40X30MM W/ LIGHT CHANNEL, RED
CT-0147	FLAT BLADE - STERNAL, MICRO, 25X25MM W/ LIGHT CHANNEL, FUCHSIA

PORT BLADES

CAT. NO.	DESCRIPTION
CT-0160	PORT BLADE CURVED, 30X50MM W/LIGHT CHANNEL, GOLD
CT-0161	PORT BLADE UPPER RIB (REVERSE CURVE), 30X50MM W/LIGHT CHANNEL, GOLD
CT-0162	PORT BLADE CURVED, SHALLOW 30X35MM W/ LIGHT CHANNEL, PURPLE
CT-0163	PORT BLADE UPPER RIB (REVERSE CURVE), SHALLOW 30X35MM W/ LIGHT CHANNEL, PURPLE

CONTACT US

TeDan Surgical Innovations, LLC

Ph: 713-726-0886

Fax: 713-726-0846

Toll Free: 877-726-0886

Email: sales@tedansurgical.com

12615 West Airport Blvd, Suite 200

Sugar Land, TX 77478

www.tedansurgical.com

TSI ©2014 LIT-0053/REVO



*For Additional Information
Contact PEAK MEDICAL
info@peakmdi.com
630-917-8234*