

VASCULAR INSTRUMENTS General Vascular Set

Call Bob
PEAK MEDICAL
630-917-8234
Bob@peakmdi.com



Clamps

AL2560.1	Gregory Profunda Clamp - Angled and Curved DeBakey Atraumatic jaws, Straight Shanks, Stainless Steel, 4.25" (11cm)
AL2562.1	Gregory Profunda Clamp - Angled and Curved DeBakey Atraumatic jaws, Straight Shanks, Stainless Steel, 5.25" (14cm)
AL2564.1	Gregory Profunda Clamp - Angled and Curved DeBakey Atraumatic jaws, Straight Shanks, Stainless Steel, 7" (18cm)
AL2566.1	Gregory Profunda Clamp - Angled and Curved DeBakey Atraumatic jaws, Straight Shanks, Stainless Steel, 7" (18cm)
AL2398.1	Wexler Multi-Purpose Clamp - 60° Angled 60mm DeBakey Atraumatic jaws, Slightly Angled Shanks, Stainless Steel, 7.25" (18.5cm)



Forceps

FL0820.2	Gerald DeBakey Forceps - Straight 1mm jaws, Flat Handle, Titanium, 7.25" (18.5cm)
FL0791.2	DeBakey Tissue Forceps - Straight 1.5mm jaws, Flat Serrated Handle, Titanium, 7.75" (19.5cm)



Needle Holders

NL0801.14	Mayo-Hegar Needle Holder - Straight jaws w/TC inserts, Stainless Steel, 7" (18cm)
NL0807.14	Ryder Needle Holder - Straight jaws w/TC inserts, Stainless Steel, 7" (18cm)
NL0723.15	Gastroviejo Micro Needle Holder - Straight jaws w/TC inserts, w/Lock, Flat Handle, Stainless Steel, 7.25" (18cm)



Scissors

SL5011.1S-0	Optimus Metzenbaum Scissors - Curved Razor edge Blades w/Blunt tips, Stainless Steel, 7" (18cm)
SL5122.1S	Diethrich-Potts Scissors - 45° Angled Razor edge Blades, Stainless Steel, 7" (18cm)
SL5125.1S	Diethrich-Potts Scissors - 125° Angled Razor edge Blades, Stainless Steel, 7" (18cm)




Retractors


TL0122.1	Weitlaner Retractor - Blunt 3x4 teeth, Stainless Steel, 6.25" (16.5cm)
TL0123.1	Weitlaner Retractor - Sharp 3x4 teeth, Stainless Steel, 6.25" (16.5cm)



VASCULAR INSTRUMENTS

Call Bob
PEAK MEDICAL
At 630-917-8234
Bob@peakmdi.com





VASCULAR INSTRUMENTS

Carotid Endarterectomy



VASCULAR INSTRUMENTS

Carotid Endarterectomy



Forceps

- FL0899.12 Ochsner Forceps - 90° Angled Very Fine serrated tips, Stainless Steel, 4.5" (11.5cm)
- FL0702.21-0 Optimus Ring tip Forceps - Straight 1mm TC coated rings, Round Handle w/Counter balance, Titanium, 7.25" (18.5cm)
- AL2262.12 DeBakey Cooley Forceps - 90° Angled Carmalt jaws, Straight Shanks, 7.5" (19cm)
- FL0804.2 DeBakey Tissue Forceps - Straight 2.5mm jaws, Flat Handle, Titanium, 7.75" (19.5cm)
- FL0807.2 DeBakey Tissue Forceps - Straight 1mm jaws, Flat Handle, Titanium, 7.75" (19.5cm)
- FL0771.21 Platform Forceps - Straight 0.5mm tips w/TC coated tying platform, Round Handle, Titanium, 8.25" (21cm)



Clamps

- AL2161.2 Micro Clamp - 60° Angled 23mm long DeBakey Atraumatic jaws, Straight Shanks, Titanium, 4" (10.5cm)
- AL2178.2 Micro Clamp - Angled DeBakey Atraumatic jaws, Slightly Curved Shanks, Titanium, 4" (10.5cm)
- AL2170.2 Micro Kitzmiller Clamp - Left Angled 23mm long DeBakey Atraumatic jaws, Straight Shanks, Titanium, 4.25" (11cm)
- AL2172.2 Micro Kitzmiller Clamp - Right Angled 23mm long DeBakey Atraumatic jaws, Straight Shanks, Titanium, 4.25" (11cm)
- AL2306.1 Wexler Baby Vascular Clamp - 45° Angled Cooley Atraumatic jaws, Straight Shanks, Stainless Steel, 5.25" (14.5cm)



Bulldog Clamps

- AL2202.2 Bulldog Shunt Clamp - 10° Angled jaws, 10 French, Titanium, 3" (7.7cm)
- AL2204.2 Bulldog Shunt Clamp - 25° Angled jaws, 10 French, Titanium, 3" (7.7cm)
- AL2206.2 Bulldog Shunt Clamp - 45° Angled jaws, 10 French, Titanium, 3" (7.7cm)
- AL2208.2 Bulldog Shunt Clamp - 60° Angled jaws, 10 French, Titanium, 3" (7.7cm)
- AL2210.2 Bulldog Shunt Clamp - 10° Angled jaws, 12 French, Titanium, 3.75" (9.5cm)
- AL2212.2 Bulldog Shunt Clamp - 25° Angled jaws, 12 French, Titanium, 3.75" (9.5cm)
- AL2214.2 Bulldog Shunt Clamp - 45° Angled jaws, 12 French, Titanium, 3.75" (9.5cm)
- AL2216.2 Bulldog Shunt Clamp - 60° Angled jaws, 12 French, Titanium, 3.75" (9.5cm)
- AL2218.2 Bulldog Shunt Clamp - 10° Angled jaws, 14 French, Titanium, 3.75" (9.5cm)
- AL2220.2 Bulldog Shunt Clamp - 25° Angled jaws, 14 French, Titanium, 3.75" (9.5cm)
- AL2222.2 Bulldog Shunt Clamp - 45° Angled jaws, 14 French, Titanium, 3.75" (9.5cm)
- AL2224.2 Bulldog Shunt Clamp - 60° Angled jaws, 14 French, Titanium, 3.75" (9.5cm)
- AL2117.1 Bulldog Kitzmiller Clamp - Left Angled 60mm long DeBakey Atraumatic jaws, Stainless Steel, 4.25" (11cm)
- AL2118.1 Bulldog Kitzmiller Clamp - Right Angled 60mm long DeBakey Atraumatic jaws, Stainless Steel, 4.25" (11cm)



Scissors

- SL5035.1S DeBakey Scissors - 60° Angled Razor edge Blades, Stainless Steel, 6.75" (17cm)
- SL5011.1ST Metzenbaum Scissors - Curved Razor edge Blades w/Tenotomy tips, Stainless Steel, 7" (18cm)
- SL5011.1S Metzenbaum Scissors - Curved Razor edge Blades w/Blunt tips, Stainless Steel, 7" (18cm)
- SL5123.1S Diethrich-Potts Scissors - 60° Angled Razor edge Blades, Stainless Steel, 7" (18cm)
- SL5129.1S Lillehei-Potts Scissors - Curved Razor edge Blades w/Blunt tips, Stainless Steel, 7" (18cm)



Needle Holders

- NL0900.21 Richards Micro Needle Holder - Straight Fine TC coated jaws, Titanium, 6" (15cm)
- NL0832.21 Ultra-Micro Needle Holder - Straight 1mm TC coated jaws, Titanium, 7.25" (18.5cm)
- NL0702.23K Castroviejo Micro Needle Holder - Straight Short TC coated jaws, w/Lock, Flat Handle, Titanium, 7.25" (18.5cm)
- NL0704.23K Castroviejo Micro Needle Holder - Straight Short TC coated jaws, w/Lock, Flat Handle, Titanium, 8.25" (21cm)



Miscellaneous

- AL1110.1 Endarterectomy Loops - 6 Gauge set, Round Handle, Stainless Steel, 5" (13cm)
- AL1058.1 Endarterectomy Spatula - 1.5mm tip, Round Handle, Stainless Steel, 7" (18cm)
- AL1050.1 Ochsner Double-Ended Dissector - 1mm/2mm tips, Round Handle, Stainless Steel, 7" (18cm)
- AL1052.1 Ochsner Double-Ended Dissector - 3mm/4mm tips, Round Handle, Stainless Steel, 7" (18cm)



Trays

- DL2218 Wexler Metal Sterilization Tray - Double level w/Pin mat and retainers, 10"x15"x3" (25x38x6.6cm)
- DL2125 Instrument Sterilization Tray - Includes Base/Mat/Lid, 15"x10"x1" (Holds 20+ Instruments)

Call Bob
PEAK MEDICAL
630-917-8234 or Bob@peakmdi.com

