

Small Patient Big Exposure



- Anatomically shaped curved rack
- Low profile hinged arms and pivoting blades
- Interchangeable blades available for pediatric and child sizes



Genesee
BIOMEDICAL, INC.

Design beyond standard.

Child, Pediatric and Neonate Sternal Retractors with Companion Atrial & Aortic Rakes

Multi-purpose Retractors for Valve or Cardiac Repair and Correction of Congenital Cardiac Abnormalities



Child Retractor with Rakes



Neonate Retractor with Rakes

Innovative Design

- Anatomic curved shape and pivoting blades distribute opening pressure differentially on the sternum.
- Curved rack with hinged arms hug patient's body unobtrusively.

Versatile Exposure

- Interchangeable rakes combined with ball clamps and H bar mounting rods produce excellent tissue exposure.
- Unique shapes of ball and pillar clamps prevent suture snagging.

Ordering Information

CATALOG NUMBER	PRODUCT DESCRIPTION
CSR-AR	Child Sternal Retractor (Curved Rack) with Atrial & Aortic Rakes
CSR-1	Child Sternal Retractor Only (Curved Rack)
PSRC-AR	Pediatric Sternal Retractor (Curved Rack) with Atrial & Aortic Rakes
PSRC-1	Pediatric Sternal Retractor Only (Curved Rack)
NSR-AR	Neonate Sternal Retractor (Curved Rack) with Atrial & Aortic Rakes
NSR-1	Neonate Sternal Retractor Only (Curved Rack)

Genesee BioMedical reserves the right to change or modify specifications without notice. U.S. Patent #5,772,583 and #6,206,828.
© 2011 Genesee BioMedical, Inc. All Rights Reserved. Printed in USA