



## CARDIAC SUMP TIPS & SUCTION TUBES

### CARDIAC SUMP TIPS



CL0802.1



CL0804.1



CL0806.1

- CL0802.1 Infant Sump Tip - 0.5" Basket, Stainless Steel, 1" (2.5cm)
- CL0804.1 Pediatric Sump Tip - 1" Basket, Stainless Steel, 1.5" (3.5cm)
- CL0806.1 Adult Sump Tip - 1.75" Basket, Stainless Steel, 2" (5cm)

### CARDIAC SUMP TIPS



CL0850.1



CL0852.1

- CL0850.1 Infant Sump Tip - 3.5mm tip, Stainless Steel, 1.4" (3.6cm)
- CL0852.1 Infant Sump Tip - Short, 3.5mm tip, Stainless Steel, 1" (2.5cm)

### CARDIAC SUMP SUCTION TUBES



CL0834.1



CL0832.1



CL0840.1

- CL0834.1 Infant 90° Sump Suction Tube - 3.5mm tip, Stainless Steel, 3.25" (8.2cm)
- CL0832.1 Infant 45° Sump Suction Tube - 3.5mm tip, Stainless Steel, 3.25" (8.2cm)
- CL0840.1 Infant 15° Sump Suction Tube - 3.5mm tip, Stainless Steel, 3.25" (8.2cm)

